



Product Datasheet

CD40 Ligand antibody (orb129552)

Description	Rabbit polyclonal to CD40 Ligand.
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CD40L(NP_000065.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	25-30KDa/17KDa
Uniprot ID	P29965
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of small intestine...



Immunofluorescence analysis of U20S cell...