



Product Datasheet

HSPB1 (phospho-S82) antibody (orb129636)


Description	Rabbit polyclonal antibody to HSPB1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	A synthetic phosphorylated peptide around S82 of human HSP27/HSPB1(NP_001531.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	28kDa
Uniprot ID	P04792
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeLa cells.



Immunofluorescence analysis of U2OS cells.



Immunoprecipitation analysis of 200ug extracts.