



Product Datasheet

KRT19 antibody (orb129657)


Description	Rabbit polyclonal antibody to KRT19
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Cytokeratin 19 (Cytokeratin 19 (KRT19))(NP_002267.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40KDa
Uniprot ID	P08727
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from m...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...