



## Product Datasheet

### NUDT6 antibody (orb129707)

|                            |                                                         |
|----------------------------|---------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to NUDT6.                             |
| <b>Reactivity</b>          | Human                                                   |
| <b>Conjugation</b>         | Unconjugated                                            |
| <b>Tested Applications</b> | ICC, IF, WB                                             |
| <b>Immunogen</b>           | Recombinant fusion protein of human NUDT6(NP_009014.2). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.        |
| <b>Concentration</b>       | 1 mg/ml                                                 |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.             |
| <b>Note</b>                | For research use only                                   |
| <b>Isotype</b>             | IgG                                                     |
| <b>Clonality</b>           | Polyclonal                                              |
| <b>MW</b>                  | 47kDa                                                   |
| <b>Uniprot ID</b>          | <b>P53370</b>                                           |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:20-50        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                         |



Western blot analysis of extracts of Jurkat and HepG2 cells.



Immunofluorescence analysis of HeLa cells.