



Product Datasheet

p53 (phospho-S9) antibody (orb129715)

Description	Rabbit polyclonal antibody to p53 (phospho-S9)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphorylated peptide around S9 of human p53(NP_000537.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	55kDa
Uniprot ID	P04637
Dilution Range	WB: 1:500-2000, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of 293...



Immunofluorescence analysis of U2OS cell...