

Product Datasheet

PI3 Kinase p110 alpha antibody (orb129718)

Description Rabbit polyclonal to PI3 Kinase p110 alpha.

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, WB

Immunogen Recombinant fusion protein of human PIK3CA(NP_006209.2).

Preservatives PBS with 0.05% proclin300,50% glycerol,pH7.3.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 110KDa


Uniprot ID **P42336**

Dilution Range WB: 1:500-2000, IHC-P: 1:50-100


Expiration Date 12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...