


Product Datasheet

PAX3 antibody (orb12972)

Description	Goat polyclonal antibody to PAX3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PAX3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 55kDa band observed in Human Brain (Frontal Cortex) and Mouse and Rat Brain lysates (calculated MW of 53.0kDa according to Human NP_852122.1 and 52.9kDa according to Mouse NP_032807.3). Recommended concentration: 1-3µg/ml.
Protein Sequence	STVHQSTIPSNPD
Clonality	Polyclonal
MW	53.0; 24.1; 22.7; 53.5; 56.0; 44.8; 45.2; 53.3
Entrez	5077
Dilution Range	ELISA: 1:64000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Brain lys...