



Product Datasheet

RAF1 (phospho-Ser621) antibody (orb129731)


Description	Rabbit polyclonal to RAF1 (phospho-Ser621).
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphorylated peptide around S621 of human Phospho-Raf1-S621(NP_002871.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	73kDa
Uniprot ID	P04049
Dilution Range	WB: 1:500-2000, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Immunoprecipitation analysis of 200ug ex...



Western blot analysis of extracts of HeL...