



Product Datasheet

RARA antibody (orb129732)


Description	Rabbit polyclonal antibody to RARA
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human RARα(NP_001019980.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60KDa
Uniprot ID	P10276
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 cell line...



IHC-P image of human esophageal cancer t...



IHC-P staining of mouse kidney tissue us...