


Product Datasheet

alpha 1 Antitrypsin antibody (orb129737)

Description	Rabbit polyclonal antibody to alpha 1 Antitrypsin
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Alpha-1 Antitrypsin (Alpha-1 Antitrypsin (SERPINA1))(NP_000286.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51KDa
Uniprot ID	P01009
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of SW480