



Product Datasheet

STAR antibody (orb129747)

Description	Rabbit polyclonal antibody to STAR
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human StAR(NP_000340.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28KDa
Uniprot ID	P49675
Dilution Range	WB: 1:500-2000, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of SK-BR-3 cell li...



Immunohistochemical analysis of formalin...