



Product Datasheet

CTDSP1 antibody (orb12976)

Description	Goat polyclonal antibody to CTDSP1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CTDSP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 30kDa band observed in Mouse Testis lysates (calculated MW of 29.3kDa according to Mouse NP_694728.1). Recommended concentration: 0.1-0.3µg/ml.. IHC: In paraffin embedded Human Tonsil shows nucleolar staining in select cells outside the follicles. Recommended concentration, 5-10µg/ml.
Protein Sequence	EAKAQDSDKIC
Clonality	Polyclonal
MW	29.2; 29.1; 29.1
Entrez	58190
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml, IHC-P: 5-10µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Testis ly...



Immunohistochemical staining of Human To...