


Product Datasheet

POU3F1 antibody (orb12982)

Description	Goat polyclonal antibody to POU3F1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	OCT6
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band observed in Human Ovary, Placenta and Umbilical Cord lysates (calculated MW of 45.5kDa according to NP_002690.3). Recommended concentration: 1-3µg/ml.
Protein Sequence	HHALHEDGHEAQL
Clonality	Polyclonal
MW	45.5
Entrez	5453
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Umbilical...