


Product Datasheet

CACNB2 antibody (orb12984)

Description	Goat polyclonal to CACNB2.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CACNB2 (aa565-579)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize all reported isoforms NP_000715.2;NP_963890.2;NP_963884.2;NP_963891.1;NP_963887.2;NP_963865.2;NP_963864.1;NP_963866.2;NP_001161417.1.
Protein Sequence	ECNKQRSRHKSKDRY
Clonality	Polyclonal
MW	68.1; 73.6; 68.2; 71.1; 69.4; 70.8; 68.8; 68.3; 67
Entrez	783
Dilution Range	ELISA: 1:4000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate u...