


Product Datasheet

Kcnc3 antibody (orb12986)

Description	Goat polyclonal antibody to Kcnc3
Species/Host	Goat
Reactivity	Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Kcnc3 / Kv3.3 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 80kDa band observed in Mouse fetal Brain lysates (calculated MW of 81.9kDa according to NP_032448.2). Recommended concentration: 0.2-0.6µg/ml. Some minor background is detected and is blocked by the immunizing peptide. We call for caution when used
Protein Sequence	QEEVIETNRADPR
Clonality	Polyclonal
MW	81.9
Dilution Range	ELISA: 1:64000, WB: 0.2-0.6 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse fetal Bra...