


Product Datasheet

KLF13 antibody (orb12992)

Description	Goat polyclonal to KLF13.
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	KLF13 (aa45-58)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 30+28kDa bands observed in Human Heart lysates (calculated MW of 31.2kDa according to NP_057079.2). Recommended concentration: 0.1-0.3µg/ml. The 28kDa is explained by partly proteolysis. Both bands were blocked by the immunizing peptide.
Protein Sequence	TLPRVEERRDGKDS
Clonality	Polyclonal
MW	31.2
Entrez	51621
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...