


Product Datasheet

TERT antibody (orb12996)

Description	Goat polyclonal to TERT.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	TERT (aa597-611)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms NP_937983.2;NP_001180305.1.
Protein Sequence	QLRELSEAEVRQHRE
Clonality	Polyclonal
MW	127; 120
Entrez	7015
Dilution Range	ELISA: 1:32000, WB: 1 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Skeletal ...