


Product Datasheet

SLC6A3 antibody (orb12998)

Description	Goat polyclonal to SLC6A3.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SLC6A3 / DAT
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 85kDa band observed in Human Brain (Substantia nigra) lysates (calculated MW of 68.5kDa according to NP_001035.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Mazei-Robison et al, J Neurosc
Protein Sequence	QLTSSTLTNPRQSP
Clonality	Polyclonal
MW	68.5
Entrez	6531
Dilution Range	ELISA: 1:4000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Substantia...