

Product Datasheet

Krox-20 (EGR2) peptide (orb13005)

Description	Krox-20 (EGR2) peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Krox-20 (EGR2)
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb381897 (Krox-20 (EGR2) antibody)
Expiration Date	6 months from date of receipt.