

Product Datasheet

Collagen 4 alpha 1 peptide (orb13012)

Description	Collagen 4 alpha 1 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Collagen 4 alpha 1
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382133 (Collagen 4 alpha 1 antibody)
Expiration Date	6 months from date of receipt.