



Product Datasheet

TSPYL2 antibody (orb13025)

Description	Goat polyclonal to TSPYL2.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	TSPYL2 (aa441-454)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 75kDa band observed in lysates of cell lines Jurkat and MOLT4 (calculated MW of 79.4kDa according to NP_071400.1). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Pancreas. Recommended concentration: 5µg/ml.
Protein Sequence	SDIDETIHDIKISD
Clonality	Polyclonal
MW	79.4
Entrez	64061
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate u...



Immunohistochemical staining of Human Pa...