



Product Datasheet

SRD5A1 antibody (orb13030)


Description	Goat polyclonal to SRD5A1.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	SRD5A1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Prostate, Liver and Small Intestine. Recommended concentration: 3.75ug/ml.
Protein Sequence	RNLRKPGDTGYK
Clonality	Polyclonal
MW	29.5
Entrez	6715
Dilution Range	ELISA: 1:1000, IHC-P: 3.75ug/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Prostate



Immunohistochemical staining of Human Small Intestine



Immunohistochemical staining of Human Liver