

Product Datasheet

C3 antibody (HRP) (orb13073)

Description	Goat polyclonal to C3 antibody conjugated to HRP
Species/Host	Goat
Conjugation	HRP
Tested Applications	ELISA
Immunogen	human Complement C3 purified from human plasma
Target	C3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 200 µg (in 200 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	187kDa
Uniprot ID	P01024
Expiration Date	12 months from date of receipt.