


Product Datasheet

GPA33 antibody (orb1311069)

Description	Mouse monoclonal antibody to GPA33
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human GPA33 (NP_005805) produced in E.coli.
Target	GPA33
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.6 kDa
Uniprot ID	Q99795
NCBI	NM_005814
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...