


Product Datasheet

PD-L1 (CD274) antibody (orb1311119)

Description	Mouse monoclonal antibody to PD-L1 (CD274)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 19-290 of human CD274 (PD-L1) (NP_054862) produced in E.coli.
Target	CD274
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	33.28 kDa
Uniprot ID	Q9NZQ7
NCBI	NM_014143
Dilution Range	WB 1:500, IHC 1:300
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...