

Product Datasheet

CD2 antibody (orb1311218)

Description	Mouse monoclonal antibody to CD2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human CD2 (NP_001758) produced in HEK293T cell.
Target	CD2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	39.3 kDa
Uniprot ID	P06729
NCBI	NM_001767
Dilution Range	IHC 1:100
Expiration Date	12 months from date of receipt.