


Product Datasheet

Tau (MAPT) antibody (orb1311322)

Description	Mouse monoclonal antibody to Tau (MAPT)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 623-758 of human MAPT(NP_058519) produced in E.coli.
Target	MAPT
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	78.7 kDa
Uniprot ID	P10636
NCBI	NM_016835
Dilution Range	IHC 1:300
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...