


Product Datasheet

LGR5 antibody (orb1311376)

Description	Mouse monoclonal antibody to LGR5
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, FC, IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 22-543 of human LGR5 (NP_003658) produced in HEK293 cell.
Target	LGR5
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	99.8 kDa
Uniprot ID	O75473
NCBI	NM_003667
Dilution Range	WB 1:200~2000, IHC 1:100~150 (recommended antigen retrieval buffer: TEE buffer pH 9.0), IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH-3T3 (Left l...