


Product Datasheet

Cytokeratin 19 (KRT19) antibody (orb1311425)

Description	Mouse monoclonal antibody to Cytokeratin 19 (KRT19)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 240-390 of human KRT19 (NP_002267) produced in E.coli.
Target	KRT19
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	43.9 kDa
Uniprot ID	P08727
NCBI	NM_002276
Dilution Range	IHC 1:100~200
Expiration Date	12 months from date of receipt.



Western blot of cell lysates (35ug) from...