


Product Datasheet

ABAT antibody (orb1311438)

Description	Mouse monoclonal antibody to ABAT
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 29-323 of human ABAT(NP_065737) produced in E.coli.
Target	ABAT
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	53.2 kDa
Uniprot ID	P80404
NCBI	NM_020686
Dilution Range	IHC 1:100~400
Expiration Date	12 months from date of receipt.



Western blot of cell lysates (35ug) from...