


Product Datasheet

ABAT antibody (orb1311441)

Description	Mouse monoclonal antibody to ABAT
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 29-323 of human ABAT(NP_065737) produced in E.coli.
Target	ABAT
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	53.2 kDa
Uniprot ID	P80404
NCBI	NM_020686
Dilution Range	IHC 1:100~200
Expiration Date	12 months from date of receipt.



Western blot of cell lysates (35ug) from...