


Product Datasheet

MRPS7 antibody (orb1311473)

Description	Mouse monoclonal antibody to MRPS7
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 38-242 of human MRPS7 (NP_057055) produced in E.coli.
Target	MRPS7
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q9Y2R9
NCBI	NM_015971
Dilution Range	WB 1:2000, IHC 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HCT116 cell lys...