


Product Datasheet

CD146 (MCAM) antibody (orb1311474)

Description	Mouse monoclonal antibody to CD146 (MCAM)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human MCAM (NP_006491) produced in HEK293T cell.
Target	MCAM
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P43121
NCBI	NM_006500
Dilution Range	WB 1:2000, IHC 1:200
Expiration Date	12 months from date of receipt.



Western blot
of cell
lysates
(35ug)
from...