


Product Datasheet

Ki67 (MKI67) antibody (orb1311511)

Description	Mouse monoclonal antibody to Ki67 (MKI67)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1160-1493 of human MKI67 (NP_002408) produced in E.coli.
Target	MKI67
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P46013
NCBI	NM_002417
Dilution Range	IHC: 1:300
Expiration Date	12 months from date of receipt.



Western Blot analysis of HT-29 and NCI-H...