


Product Datasheet

FOXP1 antibody (orb1311534)

Description	Mouse monoclonal antibody to FOXP1
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IHC, WB
Immunogen	Full length human recombinant protein of human FOXP1 (NP_116071) produced in SF9 cell.
Target	FOXP1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	75.1 kDa
Uniprot ID	Q9H334
NCBI	NM_032682
Dilution Range	IHC 1:100, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...