


Product Datasheet

AMACR antibody (orb1311542)

Description	Mouse monoclonal antibody to AMACR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Full length human recombinant protein of human AMACR (NP_055139) produced in HEK293T cell.
Target	AMACR
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	42.2 kDa
Uniprot ID	Q9UHK6
NCBI	NM_014324
Dilution Range	IHC 1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...