


## Product Datasheet

### PDSS2 antibody (orb1311567)

|                            |                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to PDSS2                                                        |
| <b>Species/Host</b>        | Mouse                                                                                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                              |
| <b>Tested Applications</b> | ChIP, IHC, WB                                                                             |
| <b>Immunogen</b>           | Full length human recombinant protein of human PDSS2(NP_065114) produced in HEK293T cell. |
| <b>Target</b>              | PDSS2                                                                                     |
| <b>Note</b>                | For research use only                                                                     |
| <b>Isotype</b>             | IgG1                                                                                      |
| <b>Clonality</b>           | Monoclonal                                                                                |
| <b>MW</b>                  | 43.9 kDa                                                                                  |
| <b>Uniprot ID</b>          | <a href="#">Q86YH6</a>                                                                    |
| <b>NCBI</b>                | <a href="#">NM_020381</a>                                                                 |
| <b>Dilution Range</b>      | IHC 1:100~200                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                           |



Western blot of cell lysates (35ug) from...