


Product Datasheet

NY-ESO-1 (CTAG1B) antibody (orb1311611)

Description	Mouse monoclonal antibody to NY-ESO-1 (CTAG1B)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human CTAG1B (NP_001318) produced in HEK293T cell.
Target	CTAG1B
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	17.8 kDa
Uniprot ID	P78358
NCBI	NM_001327
Dilution Range	IHC 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...