


Product Datasheet

CD1C antibody (orb1311627)

Description	Mouse monoclonal antibody to CD1C
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, IHC, WB
Immunogen	Full length human recombinant protein of human CD1C(NP_001756) produced in HEK293T cell.
Target	CD1C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.7 kDa
Uniprot ID	P29017
NCBI	NM_001765
Dilution Range	IHC 1:100
Expiration Date	12 months from date of receipt.



Western Blot analysis of MOLT-4 cell lys...