


Product Datasheet

CD31 (PECAM1) antibody (orb1311640)

Description	Mouse monoclonal antibody to CD31 (PECAM1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human PECAM1(NP_000433) produced in HEK293T cell.
Target	PECAM1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	82.4 kDa
Uniprot ID	P16284
NCBI	NM_000442
Dilution Range	WB 1:500, IHC 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...