

Product Datasheet

Cpt1a peptide (orb13117)

Description	Cpt1a peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Cpt1a
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb385457 (Cpt1a antibody)
Expiration Date	6 months from date of receipt.