

## Product Datasheet

### ZNF447 (ZSCAN18) antibody (orb1311724)

**Description** Mouse monoclonal antibody to ZNF447 (ZSCAN18)

**Species/Host** Mouse

**Reactivity** Human

**Conjugation** Unconjugated

**Tested Applications** ChIP, IHC, WB

**Immunogen** Full length human recombinant protein of human ZSCAN18(NP\_076415) produced in HEK293T cell.

**Target** ZSCAN18

**Note** For research use only

**Application notes** For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

**Isotype** IgG1

**Clonality** Monoclonal


**MW** 54.6 kDa

**Uniprot ID** [Q8TBC5](#)

**NCBI** [NM\\_023926](#)

**Dilution Range** IHC 1:100~200

**Expiration Date** 12 months from date of receipt.



Western blot of cell lysates (35ug) from...