


Product Datasheet

PBX1 antibody (orb1311733)

Description	Mouse monoclonal antibody to PBX1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC
Immunogen	Full length human recombinant protein of human PBX1(NP_002576) produced in HEK293 cell.
Target	PBX1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	46.4 kDa
Uniprot ID	P40424
NCBI	NM_002585
Dilution Range	IHC 1:100~200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...