


Product Datasheet

CD64 (FCGR1A) antibody (orb1311759)

Description	Mouse monoclonal antibody to CD64 (FCGR1A)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, FC, IHC, WB
Immunogen	Full length human recombinant protein of human FCGR1A(NP_000557) produced in HEK293T cell.
Target	FCGR1A
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	42.5 kDa
Uniprot ID	P12314
NCBI	NM_000566
Dilution Range	IHC 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...