


Product Datasheet

FGF21 antibody (orb1311764)

Description	Mouse monoclonal antibody to FGF21
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC
Immunogen	Human recombinant protein fragment corresponding to amino acids 29-209 of human FGF21 (NP_061986) produced in E.coli.
Target	FGF21
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	19.4 kDa
Uniprot ID	Q9NSA1
NCBI	NM_019113
Dilution Range	IHC 1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...