


## Product Datasheet

### FGF21 antibody (orb1311765)

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to FGF21                                                                                    |
| <b>Species/Host</b>        | Mouse                                                                                                                 |
| <b>Reactivity</b>          | Human                                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                          |
| <b>Tested Applications</b> | ChIP, IF, IHC                                                                                                         |
| <b>Immunogen</b>           | Human recombinant protein fragment corresponding to amino acids 29-209 of human FGF21 (NP_061986) produced in E.coli. |
| <b>Target</b>              | FGF21                                                                                                                 |
| <b>Note</b>                | For research use only                                                                                                 |
| <b>Isotype</b>             | IgG1                                                                                                                  |
| <b>Clonality</b>           | Monoclonal                                                                                                            |
| <b>MW</b>                  | 19.4 kDa                                                                                                              |
| <b>Uniprot ID</b>          | <b>Q9NSA1</b>                                                                                                         |
| <b>NCBI</b>                | <b>NM_019113</b>                                                                                                      |
| <b>Dilution Range</b>      | IHC 1:100                                                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                       |



Immunohistochemical staining of paraffin...