

Product Datasheet

LAT2 (SLC7A8) antibody (orb1311767)

Description	Mouse monoclonal antibody to LAT2 (SLC7A8)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, WB
Immunogen	Full length human recombinant protein of human SLC7A8(NP_036376) produced in HEK293T cell.
Target	SLC7A8
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	58.2 kDa
Uniprot ID	Q9UHI5
NCBI	NM_012244
Dilution Range	IHC 1:100
Expiration Date	12 months from date of receipt.

