

Product Datasheet

ALDH1L1 antibody (orb1311799)

Description	Mouse monoclonal antibody to ALDH1L1
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human ALDH1L1 (NP_036322) produced in HEK293T cell.
Target	ALDH1L1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	98.6 kDa
Uniprot ID	O75891
NCBI	NM_012190
Dilution Range	IHC 1:100, IF 1:100
Expiration Date	12 months from date of receipt.

