


Product Datasheet

XPF (ERCC4) antibody (orb1311832)

Description	Mouse monoclonal antibody to XPF (ERCC4)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC
Immunogen	Full length human recombinant protein of human XPF (NP_005227) produced in HEK293 cell.
Target	ERCC4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	104.3 kDa
Uniprot ID	Q92889
NCBI	NM_005236
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100, IHC 1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...