


Product Datasheet

ERCC1 antibody (orb1311855)

Description	Mouse monoclonal antibody to ERCC1
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, WB
Immunogen	Full length human recombinant protein of human ERCC1 (NP_973730) produced in HEK293T cell.
Target	ERCC1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	32.4 kDa
Uniprot ID	P07992
NCBI	NM_001983
Dilution Range	WB 1:500~1000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts
(35ug)...